


AccuBurn-R (Forensic) for CD/DVD Inspector

You have a disc that has just been imaged and examined. It has significant material on it and for discovery purposes must be shared with opposing counsel. Whatever your process today is for copying such discs, you will find discs that simply cannot be copied. Either the disc cannot be read by your copying system or the copy that is produced is obviously not correct. Forensic examiners are faced with this problem every day and until AccuBurn-R for CD/DVD Inspector there was no simple solution.



Technical Note

AccuBurn-R for CD/DVD Inspector writes discs using DAO "raw" mode to write P, Q and R-W subchannels to CDs. This is significantly different from how most disc copy software and hardware works and produces identical discs even under circumstances where discs cannot ordinarily be copied accurately. When copying DVD disc images it is recommended to use write-once media of the same type (+ or -) rather than write to rewritable discs.

This product, when used in conjunction with **CD/DVD Inspector** provides the ability to produce a copy of any disc quickly and accurately. Examples of difficult to copy discs:

- VCD / SVCD / XVCD video discs
- Karaoke discs
- Unfinalized drag-and-drop discs (write-once media)
- Discs with read errors
- DVD Video

Once such a disc has been imaged using CD/DVD Inspector, it can be reproduced from the image file by this product. In most cases the copy will have the same MD5 hash value as the original - the one common exception to this is unfinalised drag-and-drop discs on write-once media.

Requirements

AccuBurn-R for CD/DVD Inspector requires:

- CD/DVD Inspector V3.0 - this product is designed to work only with CD/DVD Inspector images.
- A CD and/or DVD writer. A high-grade DVD writer from Plextor is recommended as this will produce all types of discs. There can be limitations when using some older writers.
- Windows 2000, Windows XP or Windows Vista. AccuBurn-R for CD/DVD Inspector V3.0 does not support Windows 9x operating systems.
- A minimum of 512MB of RAM is recommended with Windows XP or Windows Vista.

This product is an add-on for CD/DVD Inspector and does not provide any real benefits above the consumer product without being used in conjunction with the CD/DVD Inspector product.

Data
Duplication Ltd

4 Station Approach, Wendover, Bucks HP22 6BN
Tel: 01296 621121 Fax: 01296 621125
e-mail: info@dataduplication.co.uk
www.dataduplication.co.uk

Other Forensic Software Products

AccessData Forensic Toolkit™

The Ideal Toolkit for the Forensics Specialist

Includes Forensic Toolkit, Password Recovery Toolkit, Registry Viewer and 50 Client DNA License. Includes the Oracle Database

AccessData Password Recovery Toolkit™ (PRTK)

Recover Passwords Quickly & Easily

The Password Recovery Toolkit gives you the ability to recover passwords from well-known applications and gain access to encrypted files and drives, including the encrypted file system (EFS) from Microsoft.

AccessData Distributed Network Attack

Putting Idle Time To Work

In the past, recoveries have been limited to the processing power of one machine. DNA uses the power of machines across the network or across the world to decrypt passwords.

Infinadyne CD/DVD Inspector software

Professional software for intensive analysis and extraction of data from CD-R, CD-RW and DVD media. Tailored for professionals in data recovery, forensics, and law enforcement.

Spector

Automatically Record Everything Someone Does on Their PC

Spector AUTOMATICALLY takes hundreds of screen snapshots every hour, very much like a surveillance camera. With Spector, you will be able to see EVERY instant message, e-mail, web-site visited and keystroke typed.

Other Products:-

- Forensic Tools
- Hard Disk Duplicators
- Floppy Disk Copiers
- CD/DVD Duplicators
- CD/DVD Printers

Media Duplication Services For:-

- Floppy Disk
- CD
- DVD
- Tape
- Blank Media Supply
- Standard or Customised Packaging
- Fulfilment Service